

Chemistry Paper 3 Notes

Grass

July 2022

Intro

By Grass and me. Creating the hyperef'd table of contents is pain.

Contents

1	Intro	1
2	General Experiment notes	3
3	Titration	4
3.1	Tables	4
3.2	Pipette	5
3.3	Burette + Pipette	5
3.4	Accuracy of Instruments	5
4	Dessert Qns	6
5	Temperature Related Experiments	7
6	Q.A. Test	8
6.1	Note / AFI:	8
6.2	Best/Modal Answers for Observations:	8
7	Experimental Planning	9
7.1	Things to have	9
7.2	Notes for Steps	9
7.3	Apparatus Drawing	10
8	Examples	12
8.1	Titration	12
8.1.1	Timed Assignment 2022: Paper 3 Practical	12
8.1.2	Iodine Antiseptic	14
8.1.3	Experiment on Volumetric Analysis 4 April 2022	17
8.2	Short Experiment	19
8.2.1	Timed Assignment 2022: Paper 3 Practical	19
8.2.2	Combined Chemistry Q.A. Test 2020	21
8.3	Experimental Planning	25
8.3.1	Timed Assignment 2022: Paper 3 Practical	25
8.3.2	Hydrogen Peroxide decomposition, 4 April 2022	26

General Experiment notes

1. Always read the full instructions for the whole experiment before starting
2. Use appropriate apparatus \Rightarrow Transfer with measuring cylinder if they ask for a certain volume of a solution (NEVER use beaker for measurement)
3. Reading the volume of liquid: If meniscus curved upwards, read the volume from the bottom of the meniscus. If it curves downwards, read it from the top of the meniscus.

Titration

3.1 Tables

Normal Titration Table:

Titration Number	1	2	3
Final Burette Reading / cm ³			
Initial Burette Reading / cm ³			
Volume of R / cm ³			
Best Titration Results (✓)			


Thermometric Titration Table

Final Burette Reading / cm ³	Initial Burette Reading / cm ³	Volume of R / cm ³	Highest Temperature / °C	Total Temperature Change / °C
0.00	0.00	0.00	30.0	0.0
5.00	0.00	5.00	32.5	+2.5
10.00	0.00	10.00	35.0	+5.0
:	:	:	:	:
F _n	0.00 I _n	F _n : 0.00	H _n	H _n : +30.0

Burette

Note / AFI:

1. Wash the filter funnel with the solution to be placed in the burette
2. Place the filter funnel on the mouth of the burette when filling it up
3. Remove the filter funnel after, and for all burette readings
4. Make sure the rubber tubing is filled with the solution. i.e. Ensure the tip of the burette contains no air
5. Middle of burette clip on the middle of the rubber tubing



6. Burette Readings in 2 d.p.!
7. Burette must be straightened
8. Can use formula to calculate concentration of acid A o/r base B:
Let 1 be acid A, 2 be base B

$$\frac{M_1 V_1}{n_1} = \frac{M_2 V_2}{n_2}$$
$$M_2 = \frac{M_1 V_1 n_2}{V_2 n_1}$$

Where M_i is the molar concentration of i , V_i the volume of i , and n_i the mole ratio
(in a reaction between the two substances)

9. T2CA

- (a) T - Table: 1m
- (b) 2 - 2 d.p.: 1m
- (c) C - Consistency: 1m
- (d) A - Accuracy: 2m

3.2 Pipette

Note / AFI:

1. Do not purposely 'tap' remaining liquid out!
2. Holding of pipette; Hold on to the top portion when inserting the pipette filler (NOT the bulb!)

3.3 Burette + Pipette

1. If solution is spilled on the sides of the conical flask, use the deionised water to wash it off to ensure maximum accuracy.
2. Read the burette and pipette readings at eye level.
3. Ensure there are no bubbles when filling up the burette and pipette

3.4 Accuracy of Instruments

Burette	$\pm 0.05 \text{ (cm}^3\text{)}$
Pipette	20.00 / 25.00 cm^3
Thermometer	$\pm 0.5 \text{ (}^\circ\text{C)}$
Electronic Balance	2 d.p. (g)
Measuring Cylinder	$\pm 0.5 \text{ (cm}^3\text{)}$
Stopwatch	Nearest s
Gas Syringe	Nearest cm^3

Dessert Qns

1. A student repeated the experiment and accidentally used a 20.0cm^3 pipette to measure solution Q into the flask in each titration.

The student thought he had used a 25.0cm^3 pipette

Describe and explain the effect that this would have on the percentage by mass of iodine in the solution calculated by the student.

Ans : Using 20.0cm^3 of Q instead of 25.0cm^3 will result in lower volume of R needed [1]. This will result in a lower concentration of I_2 calculated in Q and hence a lower percentage by mass of I_2 in Q . [1]

2. Explain why it is necessary to stir the mixture in the beaker before measuring the highest temperature reached by the mixture.

Ans : It is to ensure that all the metal carbonate X has reacted completely with the respective acids. [1]

3. State one **key** source of error in the results obtained in Experiments 1 to 3 (temperature experiment). Suggest **one** improvement you could make to the experiment to reduce this error.¹

Ans : Beaker can be insulated with poor heat conductor / cover with lid nested in beaker [2]

¹Experimental Error, NOT human error

Temperature Related Experiments

1. Change in temperature / enthalpy must have "+" and "-" signs

Q.A. Test

6.1 Note / AFI:

1. Always rmb to use test tube or Q.A. test – The question might not explicitly state this
2. Hold boiling tube at an angle and reduce the flame height
3. Open the bunsen burner air hole halfway if you have problems lighting it
4. Adding excess bench reagent \Rightarrow Pour some test sample away and continue adding the reagent (if the test tube is getting too full)
5. Collect solid samples with spatula (always use the small scoop)
6. Even if there is nothing seen/heard/smelled for a particular step, the observation should still be written, as "No visible Change".
7. When writing observations rmb to put all the information down. ¹
8. Some observations required are implied, but not explicitly stated. ²

6.2 Best/Modal Answers for Observations:

1. Nothing happens \Rightarrow No visible change
2. Add aqueous ammonia slowly and stir with a glass rod until no further changes are seen \Rightarrow ³ A light blue precipitate was formed, which dissolved in excess to form a deep blue solution

¹E.g.: What is the color and smell of the gas; write colorless and odorless if nothing specific seen or smelled

²E.g.: Place the boiling tube in the test-tube rack to allow its contents to settle.

\Rightarrow After waiting, rmb to add in the observation for this; like "Blue precipitate settled at the bottom of the test tube."

³No need write in the front like: "Upon adding aqueous ammonia / silver nitrate / whatever reagents"

Experimental Planning

7.1 Things to have

(in order which they should be presented)

1. Diagram
2. Approach: (Short summary of what the experiment is going to be)¹
3. Steps: 1.
4. 2.
5. :
6. n.
7. Conclusion (Look at the qns to know what to write here, i.e results that can be concluded after the expt)²

7.2 Notes for Steps

1. Be specific. Don't leave it up to interpretation. i.e. Write down the mass / volume of reagent to be added.

¹E.g.:

- (a) Approach: Comparing the effectiveness of the 3 metal oxide catalysts by measuring the volume of oxygen produced after a fixed time
- (b) Approach: Measure mass of contents before and after heating using electronic balance

²E.g.:


- (a) The metal oxide that takes the least amount of time to produce 50cm³ of oxygen gas is the most effective catalyst, while the one that takes the longest time is the least effective catalyst.
- (b) Use the formula $\frac{\text{mass before heating} - \text{mass after heating}}{\text{mass of hydrated crystals}} \times 100\%$ to find the percentage mass of water in copper (II) sulfate crystals

7.3 Apparatus Drawing


1.	Beaker		2.	Conical flask	
3.	Test tube		4.	Boiling tube	
5.	Evaporating dish		6.	Crucible	
7.	Bunsen burner		8.	Tripod stand (with wire gauze)	
9.	Retort stand		10.	Filter funnel	
11.	Gas jar		12.	Gas syringe	
13.	Pipette		14.	Burette	

15. Draw an experimental setup for the following investigations.

(a) Measure the melting point of a solid




(b) Determine the volume of oxygen gas released from thermal decomposition of mercury(II) oxide.




(c) Determine the speed of reaction by measuring the decrease in mass of a reaction mixture.

Textbook page 353



(d) Drying moist hydrogen gas produced from a reaction mixture.



Examples

8.1 Titration

8.1.1 Timed Assignment 2022: Paper 3 Practical

5

- 2 A dibasic acid, H_2X is commonly used as a reducing agent and in disinfectants. N614
You are going to determine the volume of the dibasic acid required to neutralise a fixed volume of sodium hydroxide solution and hence, deduce the molecular formula of the acid.

Read all the instructions below carefully before starting the experiment in Question 2.

Instructions

You are going to carry out a titration experiment.

P is 0.0940 mol/dm³ aqueous sodium hydroxide.

Q is aqueous solution of a dibasic acid, H_2X .

Burette Q

Pipette P (25 cm³)

- (a) Put Q into the burette.

Pipette 25.0 cm³ of P into a conical flask.

Titrate P with Q using the indicator, methyl orange. Stop the titration when the indicator turns orange. This is the end-point of the titration.

Record your titration results in the space provided, repeating the titration as many times as you consider necessary to achieve consistent results.

Results

Titration Number	1	2	3
Final burette reading (cm ³)	24.80	24.90	24.90
Initial burette reading (cm ³)	0.00	0.00	0.00
Volume of Q (cm ³)	24.80	24.90	24.90
Best titration results (✓)	✓	✓	✓

2 d.p. !!!

per lot's average : 24.90

① Methyl orange indicator → Add just 1 or 2 drops

② Q.H Test → Used beaker instead of test tube ~~NEED TO USE TEST TUBE~~

T - 1

③ Holding of pipette → ~~Hold on to the top portion when inserting the pipette filler (not the bulb)~~

2 - 1

C - 1

A - 0

④ Read instructions in advance ~~(whole set of instructions)~~

③

[5]

⑤ Use appropriate apparatus → Transfer with measuring cylinder if they fail to certain volume of a solution.

(NEVER use beaker for measurement)

CHAPTER 8. EXAMPLES

13

- (b) From your titration results, obtain an average volume of Q to be used in your calculations. Show clearly how you obtained this volume.

~~if you must show working~~

$$\frac{14.40 + 14.90}{2} = 14.65 \text{ cm}^3$$

average volume of Q 14.65 cm³ [1]

- (c) 1 mole of the acid H₂X reacts with 2 moles of sodium hydroxide. Using your results, calculate the concentration, in mol/dm³, of H₂X in Q.

$$\begin{aligned} \text{Moles of NaOH} &: \text{H}_2\text{X} \\ 2 &: 1 \\ 0.01235 \text{ mol} &\rightarrow 0.001175 \text{ mol} \\ \text{Conc. of H}_2\text{X in Q} &= \frac{0.001175}{(1000)} \times 1000 \text{ mol/dm}^3 \\ &= 0.001175 \text{ mol/dm}^3 (1) \cdot f \\ \text{concentration of H}_2\text{X in Q} &= 0.001175 \text{ mol/dm}^3 \quad \text{eef.} \end{aligned}$$

[2]

- (d) Assuming that the concentration, in g/dm³, of H₂X is 3.53 g/dm³, use your answer in (c) to calculate the relative molecular mass of H₂X.

$$\begin{aligned} \text{M}_r \text{ of H}_2\text{X} &= \frac{3.53}{(0.001175)} \\ &= 3000 \text{ g/mol.} \end{aligned}$$

relative molecular mass of H₂X 3000 eef. [1]

- (e) The expected molecular formula of H₂X is determined to be either H₂SO₃ or H₂SO₄. Based on your answer in (d), identify the molecular formula of H₂X.

$$\begin{aligned} \text{M}_r \text{ of H}_2\text{SO}_3 &= 2(1) + 3(16) \\ &= 74 \\ \text{M}_r \text{ of H}_2\text{SO}_4 &= 2(1) + 4(16) \\ &= 78 \end{aligned}$$


*Your answer does not contain either H₂SO₃ or H₂SO₄.
Lucky guess!!*

molecular formula of H₂X H₂SO₃ [3]

- (f) A student repeated the experiment and accidentally used a solution of sodium carbonate, Na₂CO₃, instead of sodium hydroxide in the titration. The solution of sodium carbonate is of the same concentration as P.

Would the volume of Q required be higher or lower than the titre values obtained in (b)? Explain why.

*Volume of Q required would be higher.
Number of moles of Q required would be higher as 1 mole of Na₂CO₃ reacts with 1 mole of H₂X in Q.*



N : H


(3)

8.1.2 Iodine Antiseptic

Name: Yu Shau Hung (14) Class: 4-I Date: _____ Mark: _____ / 23

- 2 Iodine is widely used as an antiseptic.

You are going to determine the concentration of iodine in an antiseptic solution by titration with aqueous sodium thiosulfate, $\text{Na}_2\text{S}_2\text{O}_3$, using starch solution as an indicator.



Read all the instructions below carefully before starting the experiment in Question 2.

Instructions

You are going to carry out a titration experiment.

Q is an antiseptic solution containing iodine.

R is 0.160 mol/dm³ aqueous sodium thiosulfate.

- (a) (i) Put **R** into the burette.

Pipette 25.0 cm³ of **Q** into a flask.

Add **R** from the burette until the red-brown colour fades to pale yellow, then add a few drops of the starch solution. This will give a dark blue solution.

Continue adding **R** slowly from the burette until one drop of **R** causes the blue colour to disappear, leaving a colourless solution.

Record your titration results in the space provided, repeating the titration as many times as you consider necessary to achieve consistent results.

Results

Titration number	1	2	3
Final burette reading /cm ³	23.20	46.50	32.50
Initial burette reading /cm ³	0.00	24.20	10.00
Volume of R /cm ³	23.20	22.30	22.50
Best titration results (✓)		✓	✓

[5]

- SYNTHETIC / ORGANIC
- (ii) From your titration results, obtain an average volume of R to be used in your calculations.
Show clearly how you obtained this volume.

$$\frac{21.10 + 22.50}{2} = 22.40$$

average volume of R 22.40 cm³

- (b) (i) R is 0.160 mol/dm³ aqueous sodium thiosulfate.

Using your results, calculate the concentration, in mol/dm³, of iodine in Q.

$$\begin{aligned}
 & \text{mol of R} && \text{Na}_2\text{S}_2\text{O}_3 : \text{I}_2 \\
 & = 0.160 \times \frac{22.40}{1000} && 2 : 1 \\
 & && 0.003584 \text{ mol} \rightarrow 0.001792 \text{ mol} \\
 & = 0.003584 \text{ mol} && \text{concentration of I}_2 \text{ in Q} \\
 & && = \frac{0.001792}{(1000)} \\
 & && = 0.0001792 \text{ mol/dm}^3 \\
 & && = 0.0001792 \text{ mol/dm}^3 \text{ (3.s.f.)} \\
 & && \text{concentration of iodine in Q} 0.0001792 \text{ (3.s.f.)} \text{ mol}
 \end{aligned}$$

- (ii) Using your answer from (i), calculate the mass of iodine in 1 dm³ of Q. [A: 1, 127]

$$\begin{aligned}
 & \text{mol of Q} && \text{mass of I}_2 \text{ in Q} \\
 & = 0.0001792 \times 1 && = 0.0001792 \times 2(127) \\
 & = 0.0001792 \text{ mol} && = 0.0001792 \text{ g} \\
 & && = 0.0001792 \text{ g} \\
 & && \text{mass of iodine in } 1 \text{ dm}^3 \text{ of Q} 0.0001792 \text{ g (3.s.f.)}
 \end{aligned}$$

- (iii) Given that 1 cm³ of Q has a mass of 1.00 g, calculate the percentage by mass of the solution.

$$\begin{aligned}
 & \cancel{\text{conc. of Q in g/dm}^3} && \text{mass of } 1 \text{ dm}^3 \text{ of Q} && \% \text{ by mass of I}_2 \text{ in Q} \\
 & \cancel{\frac{1}{(1000)}} && = 1 \times 1000 && = \frac{0.0001792}{1000} \times 100\% \\
 & = 1000 \text{ g/dm}^3 && = 1000 \text{ g} && = 0.0001792\% \\
 & && && = 0.0001792 \text{ g} \\
 & && && = 0.0001792 \text{ g (3.s.f.)}
 \end{aligned}$$

percentage by mass of iodine in Q 0.0001792% (3.s.f.)

- (c) A student repeated the experiment and accidentally used a 20.0 cm^3 pipette to measure solution Q into the flask in each titration.

The student thought he had used a 25.0 cm^3 pipette.

Describe and explain the effect this would have on the percentage by mass of iodine in the solution calculated by the student.

The student would calculate a lower concentration of iodine in Q as the mol of R, and hence Q, would decrease proportionately when using a 20.0 cm^3 pipette instead, but the student will calculate it using a higher volume of 25.0 cm^3 instead of the correct volume of 20.0 cm^3 . So the percentage of iodine calculated, and hence the percentage by mass of iodine, will decrease. [2]


A: Using 20.0 cm^3 of Q instead of 25.0 cm^3 will result in lower volume of R [Total: 14 marks] [1]

This will result in lower concentration of I_2 calculated in Q, and hence a lower percentage by mass of I_2 in Q. [1]

8.1.3 Experiment on Volumetric Analysis 4 April 2022

Redox Titration

Aqueous hydrogen peroxide, H_2O_2 , can be prepared in the laboratory by adding dilute sulfuric acid to barium peroxide, BaO_2 , and filtering the resulting mixture to remove insoluble barium sulfate.



$\rightarrow P$: contains excess $H_2SO_4(aq)$.

Conc. of P ?

V of $P = 10.0\text{ dm}^3$

P is an acidified solution of hydrogen peroxide prepared in this way. You are to determine its concentration and the mass of barium peroxide required to prepare 10.0 dm^3 of P .

The concentration of hydrogen peroxide can be determined by titrating it with potassium manganate(VII). No indicator is necessary since the products of the reaction are almost colourless and one drop of potassium manganate(VII) in excess produces an easily seen pale pink colour.

Manganate peroxide?

(conc.?)

Q is 0.0200 mol/dm^3 potassium manganate(VII), $KMnO_4$. ($P \xrightarrow{\text{Pip}} P \Rightarrow$ Acidified H_2O_2 (conc.?)

(Bulb) $Q \Rightarrow 0.02\text{ mol/dm}^3 KMnO_4$ at 25.0 cm^3

- (a) Fill the burette with solution Q . You may find it easier to read the top of the meniscus as the colour of Q is so intense.

Pipette a 25.0 cm^3 (or 20.0 cm^3) portion of P into a flask and titrate with Q . At first, the purple colour disappears rapidly. As the titration proceeds, this disappearance is less rapid. At the end-point, one drop of Q produces a pink colour that does not disappear on swirling.

colourless \Rightarrow Purple pink

Record your titration results in the space provided, repeating the titration as many times as you consider necessary to achieve consistent results.

Results

Titration No.	1	2	3
Initial burette reading / cm^3	0.00	11.20	16.80
Final burette reading / cm^3	26.40	37.20	43.00
Volume of Q / cm^3	26.40	26.00	26.10 26.20
Best titration result (✓)		✓	✓

T - 0
Z - 1
C - 1
A - 0

- (b) From your titration results, obtain an average volume of Q to be used in your calculations.

Show clearly how you obtained this volume.

$$\frac{26.08 + 26.16}{2} = 26.05 \text{ cm}^3$$

ref.
average volume of Q ... 26.05 cm³

- (c) Q is 0.0200 mol/dm³ potassium manganate(VII).

Five moles of hydrogen peroxide react with two moles of potassium manganate(VII)

Using your results from (b), calculate the concentration, in mol/dm³, of the hydrogen peroxide in P.

$$\begin{aligned} \text{mol of Q} &= \frac{26.05}{1000} \times 0.0200 \\ &= 0.000521 \text{ mol} \end{aligned}$$

$\text{H}_2\text{O}_2 : \text{KMnO}_4$
5 : 2
 $\frac{0.000521 \text{ mol}}{2} \leftarrow 0.0002605 \text{ mol}$

$\text{Conc. of H}_2\text{O}_2 \text{ in P}$
 $= \frac{0.0002605}{\left(\frac{26}{1000}\right)}$
 $= 0.0521 \text{ mol/dm}^3$

$$\begin{aligned} I &= \text{H}_2\text{O}_2, \quad 2 = \text{KMnO}_4 \\ \frac{M_1 V_1}{n_1} &= \frac{M_2 V_2}{n_2} \\ \frac{M_1 \times 25}{5} &= \frac{0.02 \times 27.50}{2} \\ M_1 &= 0.0521 \text{ mol/dm}^3 \end{aligned}$$

Concentration of hydrogen peroxide in P is 0.0521 mol/dm³.

[2]

- (d) Using your answer in (c), calculate the mass of barium peroxide needed to prepare 10.0 dm³ of solution P.
[Ar: O, 16; Ba, 137]

$$\begin{aligned} \text{mol of H}_2\text{O}_2 \text{ in } 10.0 \text{ dm}^3 \text{ of P} &= 0.0521 \times 10.0 \\ &= 0.521 \text{ mol} \end{aligned}$$

$\text{BaO}_2 : \text{H}_2\text{O}_2$
1 : 1
 $0.521 \text{ mol} \leftarrow 0.521 \text{ mol}$

$$\begin{aligned} \text{Mass of BaO}_2 &= 0.521 \times [137 + 2(16)] \\ &= 88.049 \text{ g} \end{aligned}$$

(88.049)

Mass of barium peroxide required is 88.0 g.

[2]

3ef.

8.2 Short Experiment

8.2.1 Timed Assignment 2022: Paper 3 Practical

1 The reaction between acid and metal carbonate is exothermic. The reaction can be represented by the equation:

$$2\text{H}^+(\text{aq}) + \text{CO}_3^{2-}(\text{aq}) \rightarrow \text{CO}_2(\text{g}) + \text{H}_2\text{O}(\text{l})$$

You are going to investigate the relative strength of acid solutions A, B and C, when a fixed mass of metal carbonate X is dissolved in each of the acid solution in excess. All the acids have the same concentration and are monobasic acids of different acid strength.

You are also going to suggest the identity of the cation in metal carbonate X.

Read all the instructions below carefully before starting the experiment in Question 1.

Instructions

You are going to carry out four experiments.

You are provided with acid solutions A, B and C and three vials of solid metal carbonate X of the same mass. The masses are pre-weighted and labelled for you.

(a) **Experiment 1**

Using a measuring cylinder, transfer about 40 cm^3 of acid solution A into a 100 cm^3 beaker. Add a vial of solid X in the beaker, and record it in the table. Measure the highest temperature reached and record it in the table. Wash the thermometer and beaker with distilled water.

Experiment 2 and 3

Repeat the procedure in Experiment 1, using acid solutions B and C with the two remaining vials of solid X. Keep the resulting mixture from Experiment 3 for Experiment 4. This mixture will be known as solution Y.

Complete the table by calculating the change in temperature for each experiment.

Results:

experiment	acid solution	initial temperature / °C	highest temperature reached / °C	change in temperature / °C
1	A	31.0	34.0	+3.0
2	B	31.0	32.0	+1.0
3	C	31.0	35.0	+4.0

Keep resulting mixture

T - Data showing increasing trend - 1m T - 1
 R - Temperature reading within $\pm 0.5^\circ\text{C}$ of previous's results - 1m R - 1
 D - All values in 1 decimal place - 1m D - 1
 S - Positive sign indicated for change in temp - 1m S - 0

3

Experiment 4

greenish B/g
blue soln

To a portion of solution Y, add aqueous ammonia until no further change is observed. Record your observations and suggest the identity of the cation present.

When aqueous ammonia is added, the green precipitate is formed which is insoluble in excess ammonia. However, the cation present is Fe^{2+} from $\text{Fe}^{2+}(\text{aq}) + 2\text{OH}^-(\text{aq}) \rightarrow \text{Fe(OH)}_2(\text{s})$. Give precipitate dissolves in excess aqueous ammonia to form a dark blue solution. $\text{Fe}^{2+}(\text{aq}) + 2\text{OH}^-(\text{aq}) \rightarrow \text{Fe(OH)}_2(\text{s})$

(b) Give two observations made when solid X reacts with the acid solutions and explain each observation.

observation 1 Effect of strong acid, i.e. gas is evolved. (carbon dioxide)
 explanation The following diagram shows when strong acid solution react with metal carbonate X, there is a temperature change due to evolution of carbon dioxide bubbles.

observation 2 The solution turns darker
 explanation The reaction between solid X and the acid solution are exothermic.
 In, thermal energy is released due to the reaction, causing the solution to be hotter and temperature to increase.

(c) (i) Using the given formula as shown, calculate the enthalpy change, ΔH for experiments 1, 2 and 3 in the space provided.

$\Delta H (\text{kJ}) = \text{total volume of mixture } (\text{cm}^3) \times 0.0042 \times \text{temperature change } (^{\circ}\text{C})$

ΔH for experiment 1
 $= 40 \times 0.0042 \times 3.0$
 $= 0.504 \text{ kJ}$

ΔH for experiment 2
 $= 40 \times 0.0042 \times 1$
 $= 0.168 \text{ kJ}$

ΔH for experiment 3
 $= 40 \times 0.0042 \times 4$
 $= 0.172 \text{ kJ}$

-1 no sign!

(ii) Hence, arrange the relative strength of acid solutions A, B and C in ascending order, starting with the acid solution of the lowest strength. Explain your answer.

Relative strength of acid solutions is in descending order, i.e. $\text{H}_2\text{A} > \text{B} > \text{C}$.

Since the acid solutions are the first reactant, and the molar concentration is fixed, the temperature change depends on the relative strength of each acid, meaning the higher the relative strength of each acid, the higher the ΔH and the temperature.

- (d) Explain why it is necessary to stir the mixture in the beaker before measuring the highest temperature reached by the mixture.

~~It is to ensure that all the metal carbonate X has reacted with the respective acid.~~

~~(Also increases speed of reaction between X and the acids)~~

- (e) State one key source of error in the results obtained in Experiments 1 to 3. Suggest one improvement you could make to the experiment to reduce this error.

~~error thermal energy will be lost to the surroundings easily as the top of the beaker is produced exothermically~~


~~Expt error, not human error]~~

~~improvement cover the beaker quickly after the reaction has completed (and/or) wrap the beaker with some insulating materials, like a cloth.~~

- (f) If the same experiments were repeated using an excess of metal oxide instead of metal carbonate, the temperature change would be positive. Suggest why this is so.

~~Metal oxides are basic oxides. So, when they react with the acid solutions, they neutralise each other, where $\Delta H < 0$, releasing thermal energy into the surroundings, thereby creating a positive temperature change.~~

[2]



(Q)

8.2.2 Combined Chemistry Q.A. Test 2020



8

For Examiner's Use

2 You are provided with:

- two samples of solid P
- a sample of solution Q
- two samples of solution R, labelled R₁ and R₂
- a sample of solution I (screened methyl orange indicator).

Q is a salt and solutions R₁ and R₂ have different concentrations.

Carry out the following tests. You should test and identify any gases evolved. Carefully record your observations.

The volumes given below are approximate and should be estimated rather than measured unless instructed otherwise.

test	observations
<p>(a) Add about 5cm depth of limewater into a clean test-tube.</p> <p>Remove the stopper from one of the boiling tubes containing solid P.</p> <p>Heat this boiling tube containing solid P gently for about one minute.</p> <p>While you are heating the boiling tube containing solid P, withdraw several samples of gas from the boiling tube using a dropper pipette.</p> <p>Do not allow the dropper pipette to touch the boiling tube.</p> <p>Bubble each sample of gas into the limewater.</p>	<p>P: CO_3^{2-} ✓</p> <p>The samples of gas formed a white precipitate with the limewater.</p> <p>B: No visible change (during heating)</p> <p>Colourless and odourless gas evolved, produced a white precipitate in limewater.</p>
<p>(b) Add about 1cm depth of solution Q into the (other) boiling tube containing solid P. Shake the mixture.</p> <p>Place the boiling tube in the test-tube rack to allow its contents to settle.</p>	<p>P: Cu^{2+} ✓</p> <p>When Q was added to P, a light blue precipitate was formed.</p> <p>B: Blue precipitate produced Bubbling occurs Colourless and odourless gas produced a white precipitate when bubbled through limewater. Blue precipitate settled at bottom of test-tube</p>



Small
 ↓ NH
 ↓ CO₂
 ↓ O₂ or
 ↓ H₂ or
 SO₂

* 0012511521209 *

9

	test	observations
(c)	<p>Add about <u>1cm</u> depth of <u>solution Q</u> into a clean test-tube.</p> <p>To this test-tube, add about <u>1cm</u> depth of dilute nitric acid.</p> <p>To the same test-tube, add about <u>1cm</u> depth of barium nitrate solution. $Q: SO_4^{2-} X$</p>	<p>When barium nitrate was added, the solution remained blue. ✓</p> <p>OR: No visible reaction [1]</p>
(d)	<p>Add about <u>1cm</u> depth of <u>solution Q</u> into a clean test-tube.</p> <p>To this test-tube add about <u>1cm</u> depth of dilute nitric acid.</p> <p>To the same test-tube, add about <u>1cm</u> depth of silver nitrate solution. $Q: Cl^- \checkmark$</p> <p>Then add aqueous ammonia slowly and stir with a glass rod until no further changes are seen.</p> <p>$Q: Cu^{2+} \checkmark$</p>	<p>When silver nitrate was added, a white precipitate was formed.</p> <p>Upon adding aqueous ammonia, a light blue precipitate was formed, which dissolved in excess aqueous ammonia to form a deep blue solution [3]</p>
(e)	<p>Add about <u>1cm</u> depth of <u>solution Q</u> into a clean test-tube.</p> <p>To this test-tube, add aqueous sodium hydroxide slowly with shaking, until no further changes are seen.</p>	<p>Upon adding aqueous sodium hydroxide, a light blue precipitate was formed which is insoluble in excess aqueous sodium hydroxide. ✓</p> <p>[1]</p>

(f) (i) Using your observations, deduce the identity of salt Q.

For the cation and for the anion in salt Q, justify your deductions using evidence from your observations.

salt Q is ...Copper (II) chloride.....

^{(b) In (c)}

justification for the cation ...Upon adding aqueous ammonia, a light blue precipitate formed.....

use either
reaction from testing
NaOH or NH_3

.....which dissolved in excess^{(a) aqueous ammonia}.....to form a colourless solution. And when aqueous sodium hydroxide^{(b) sodium hydroxide}.....

.....was added, a light blue precipitate formed which is insoluble in excess^{(c) dilute nitric acid}.....

justification for the anion ...When silver nitrate was added^(d), a white precipitate was
.....formed, which must be Cl^-

[2]

(ii) Suggest one conclusion about P that can be drawn from your observations. Support your conclusion with evidence from your observations.

conclusion ...P contains carbonate anions.....

^{(a), (b)}

evidence ...When P was heated^(c), the gas evolved reacted with limewater to form a
.....white precipitate, meaning it was likely a carbonate salt which decomposed to give the
.....carbon dioxide gas that reacted with limewater.....

[1]

- (g) In this question, you are going to compare the concentrations of two solutions of a basic compound, R_1 and R_2 .

<p>(i)</p> <p>Using a measuring cylinder, transfer 2.0 cm^3 of solution R_1 into the test-tube labelled R_1.</p> <p>Using a different measuring cylinder, transfer 2.0 cm^3 of solution R_2 into the test-tube labelled R_2.</p> <p>Add about 6 drops of solution I (screened methyl orange indicator) to both test-tubes.</p> <p>Using a dropper pipette, add dilute nitric acid, drop by drop with shaking, to R_1 until the indicator changes from green to purple.</p> <p>Count and record the number of drops of acid added to solution R_1.</p> <p>Using the same dropper pipette, add dilute nitric acid, drop by drop with shaking, to R_2 until the indicator changes from green to purple.</p> <p>Count and record the number of drops of acid added to solution R_2.</p>	<p>number of drops of acid added to R_1 = 9</p> <p>number of drops of acid added to R_2 = 6</p>
---	---

[2]

- (ii) Use your results from (g)(i) to deduce which of the solutions R_1 or R_2 is more concentrated. Explain your answer.

The more concentrated solution of R is ... solution R_1 . ✓

explanation ... Given the same volume of R_1 and R_2 , R_1 took more acid to reach the endpoint of the titration, suggesting there was a greater number of moles of the basic compound in R_1 , and hence making R_1 more concentrated than R_2 . [1]

OR: Larger volume of nitric acid is needed to neutralise the more concentrated R_1 .

[Total: 15]

8.3 Experimental Planning

8.3.1 Timed Assignment 2022: Paper 3 Practical


3 Copper(II) sulfate crystals contain water of crystallisation which may be removed by heating.

- (a) Given a sample of copper(II) sulfate crystals and the usual apparatus found in a school laboratory, plan an experiment to find the percentage by mass of water in copper(II) sulfate crystals.

You should include the measurements you would take and show how you would use your results to calculate the percentage by mass of water in copper(II) sulfate crystals. You should describe or draw a diagram of the apparatus that may be used to remove the water.

Diagram

Measure mass of crystals before and after heating using electronic balance.



Heat to ensure that all water has been removed. Approach: Measure and record the change in mass after heating hydrated copper(II) heat to dryness/ heat to constant mass / heat until solid turns white. In state crystals till they turn anhydrous. Then calculate the percentage by mass of water in $\text{CuSO}_4 \cdot y\text{H}_2\text{O}$. Steps: 1. Set up the experiment as shown in the diagram, but with an electronic balance where the mass of the empty crucible and the crucible with 5g of CuSO_4 crystals are measured with the electronic balance. 2. Use the Bunsen burner to heat the hydrated CuSO_4 crystals till they become anhydrous. 3. Measure the mass of CuSO_4 with the anhydrous CuSO_4 left using the electronic balance. 4. Calculate change in mass of the CuSO_4 : Mass of crystals with 5g of CuSO_4 crystals - mass of crystals with only water in the copper(II) sulfate crystals. 5. Using the formula $\frac{\text{mass of CuSO}_4 \text{ (Step 4)}}{\text{mass of CuSO}_4 \text{ (Step 4)}} \times 100\%$, find the percentage by mass of water in the copper(II) sulfate crystals.

- (b) The formula for copper(II) sulfate crystals is $\text{CuSO}_4 \cdot y\text{H}_2\text{O}$ where y is the number of moles of water of crystallisation in one mole of crystals. A student does an experiment and finds that $y = 4$. The correct value of y for her sample is 5.

Suggest an error in her experiment that would result in this difference. Explain how this error would lead to a lower value of y .

You can assume that all her measurements were read and recorded correctly and that her calculation was correct.

The ^{hydrated} copper(II) sulfate crystals did not fully turn anhydrous, and some water

of crystallisation was left after heating in the crystals. Therefore, a higher

mass of supposedly anhydrous CuSO_4 crystals is measured, resulting in a lower

change in mass of the crystals, which is the mass of water in the crystals.

OR → Recalculated. So, the percentage by mass of water in the hydrated CuSO_4 crystals also decreased, resulting in a lower number of moles of water, and therefore a lower mole ratio, and y calculated.


if water

[Total: 6]

5


8.3.2 Hydrogen Peroxide decomposition, 4 April 2022

- (f) Hydrogen peroxide decomposed easily, as shown in the equation below.



[4]

The reaction may be catalysed by several different metal oxides including copper(II) oxide, iron (III) oxide and manganese (IV) oxide. Given a supply of aqueous hydrogen peroxide, the powdered metal oxides and the apparatus normally found in a school chemistry laboratory, describe an experiment which would enable you to decide which of the three metal oxides is the most effective catalyst and which is the least effective. Include a labelled diagram in your answer.



** End of Practical **

2

Approach: Effectiveness of the 3 metalⁿ catalysts can be compared by measuring the volume of oxygen gas produced after a fixed time.

1. Set up the experiment as shown above, adding 1.0g of powdered copper (II) oxide at the last.
2. Upon adding the 1.0g of powdered copper (II) oxide, immediately seal the conical flask using the stopper and start the stopwatch.
3. Stop the stopwatch when 50 cm³ of oxygen gas has been collected.
4. Repeat the experiment twice, using powdered iron (III) oxide or manganese (IV) oxide each once in place of powdered copper (II) oxide.
5. The metal oxide that takes the least amount of time to produce 50 cm³ of oxygen gas is the most effective catalyst, while the one that takes the longest time is the least effective catalyst.